

PERFORMANCE DATA SHEET

2610

Monoclonal anti-human Lambda Light Chain/FITC***mAb name/Clone:** ICO106**Isotype:** Mouse IgG1**Immunogen:** Human Ig**CATALOG#:** 138-040**QUANTITY:** 120 tests**VOLUME IN VIAL:** 0.2ml**WORKING DILUTION:** 1:50 (or use 1.6µl of concentrated stock per 5×10^5 -cell test)

INFORMATION: Human immunoglobulins are glycoproteins composed of two disulfide-bonded heavy (H) chain subunits, each of which is linked by interchain disulfide bonds to a light (L) chain forming a tetramolecular complex. There are five classes of immunoglobulins, designated IgG, IgA, IgM, IgD and IgE, which are defined by differences in the constant region of H chains. L chains are divided into kappa or lambda classifications based on structural antigenic differences. All classes of immunoglobulins have been found on the cell surface of B lymphocytes where they function as antigen receptors to elicit antigen-dependent proliferation and secretion of antigen specific soluble circulating antibodies. Antibody ICO106 recognizes cell surface lambda chain on B cells

References: 1. Basic and Clinical Immunology, Seventh edition (D. P. Sites & A. I. Terr, eds.) Appleton & Lange., Norwalk, CT (1991).

STORAGE CONDITIONS: Store at 2 - 5°C. Freeze/Thawing is not recommended. Protect from light.

PRODUCT STABILITY: Product should retain activity for at least 12 months after shipping date when stored as recommended. Ship Date: _____

BUFFER: 50 mM Sodium Phosphate pH 7.5, 100 mM Potassium Chloride, 150mM NaCl, 5% Glycerol, 0.2% BSA, 0.04% NaN₃ (as a preservative).

PRODUCTION: Protein A purified antibody from tissue culture supernatant was reacted with FITC. Unconjugated FITC was separated from antibody/FITC conjugate by desalting column. The antibody/FITC conjugate is at 0.5 mg/ml with a Fluorescein/IgG molar ratio of 7.9.

PERFORMANCE: Five $\times 10^5$ ficoll prepared human peripheral blood mononuclear cells were pre incubated 5 minutes with 20 µl of 250 µg/ml human IgG after which they were incubated 45 minutes on ice with 80 µl of anti-Lambda/FITC at a 1:50 dilution (10 µg/ml). They were then washed three times, fixed and analyzed by FACS. Cells stained positive with a mean shift of 1.23 log₁₀ fluorescent units when compared to a Mouse IgG1/FITC negative control (Catalog # 278-040) at a similar concentration. Binding was blocked when cells were pre incubated 10 minutes with 20 µl of 0.5 mg/ml anti-lambda antibody (Catalog #138-020).

* Research Use Only. Not for use in Diagnostic procedures.

Binding of anti-Lambda/FITC to human PBMC

