PERFORMANCE DATA SHEET 3241 Monoclonal anti-human CD54(ICAM-1)/FITC*



mAb name/Clone: **8.4A6** *Isotype:* Mouse IgG1 *Immunogen:* Human endothelial cells

CATALOG#: 206-040 QUANTITY: 120 tests WORKING DILUTION: 1:50 (or use 1.6ul of concentrated stock per 5 x 10⁵-cell test)

INFORMATION: Human CD54 (ICAM-1) mediates cell adhesion by binding to the integrins CD11a/CD18 (LFA-1) and to CD11b/CD18 (Mac-1). CD54 expression on resting peripheral blood leukocytes is weak but is upregulated on activated T and B lymphocytes and monocytes. Antibody 8.4A6 recognizes the D2 domain of the CD54 molecule of 90 kd. Antibody 8.4A6 inhibits CD54 binding to LFA-1. *References:* A.R. Berendt, et al, (1992) Cell **68**: 71-81. Leukocyte Typing V (S.F. Schlossman, et al, eds.) Oxford University Press, Oxford, (1995) p. 1548-1550. P.L. Reilly, et al, (1995) J Immunol **155**: 529-532.

STORAGE CONDITIONS: Store at 2 - 5^oC. Freeze/Thawing is not recommended. Protect from light.

PRODUCT STABILITY: Product should retain activity for at least 12 months after shipping date when stored as recommended. Ship Date:_____


BUFFER: 50 mM Sodium Phosphate pH 7.5, 100 mM Potassium Chloride, 150mM NaCl, 5% Glycerol, 0.2% BSA, 0.04% NaN₃ (as a preservative).


PRODUCTION: Protein A purified antibody from tissue culture supernatant was reacted with FITC. Unconjugated FITC was separated from antibody/FITC conjugate by desalting column. The antibody/FITC conjugate is at approximately **0.25 mg/ml** with a

Fluorescein/IgG molar ratio of 9.1.

PERFORMANCE: Five x 10^5 cultured **Raji** cells were washed and incubated 45 minutes on ice with 80 µl of anti-CD54/FITC at a **1:50** dilution (5 µg/ml). They were washed three times, fixed and analyzed by FACS. Cells stained positive with a mean shift of **1.52** log₁₀ fluorescent units when compared to a Mouse IgG1/FITC negative control (Catalog # 278-040). Binding was blocked when cells were pre incubated 10 minutes with 20 µl of 0.5 mg/ml anti-CD54 (8.4A6) antibody (Catalog #206-020).

* Research Use Only. Not for use in Diagnostic procedures..





Ancell Corporation P.O. Box 87 Bayport, MN 55003-0087 USA Phone: Toll free 800-374-9523 or 651-439-0835 Fax: 651-439-1940